Bipolar coagulation forceps

Symbols on product and packages

Symbol	Explanation
	CAUTION: general warning symbol CAUTION: see documentation supplied with the product.

Intended use

Bipolar forceps from Ulrich Swiss are intended for use in all areas of open surgery. They are designed to grasp, manipulate and coagulate tissue. They are connected through a suitable bipolar cable with the bipolar output of a HF generator and may be used only with bipolar coagulation current. Ulrich bipolar forceps must be operated with the following parameters: Frequency range between 300 kHz and 1,000 kHz; maximum generator operating voltage of 600Vp.

The surgeon must determine the appropriate type of treatment for each individual case in cooperation with the internist and the anaesthetist. Typical HF surgical interventions include

- Dermatology: hair and wart removal
- ENT: polyp removal, bipolar hemostasis
- Gynaecology: electro-conization of the portio, sterilization by bipolar coagulation of the fallopian tube
- Urology: transurethral resection of the prostate, bladder tumor or bladder wall papillomas; cutting and coagulation in transabdominal surgery
- Laparoscopy and cystoscopy: Removal of polyps
- General surgery: primary sections, separating tissue, plastic surgery, hemostasis, coagulation
- Neurosurgery: bipolar coagulation

Contraindications

Contraindications for the use of electrosurgical instruments are given when the physician's judgement would find their use contrary to the best interest of the patient. The Ulrich Swiss bipolar forceps have not been shown to be effective for tubal sterilization or tubal coagulation for sterilization procedures and should not be used for this purpose.

No or Limited HF surgical interventions should take place

- for patients with implantable cardiac pacemakers or pacemakers electrodes
- in areas where flammable or explosive materials are present, for example in the gastro-intestinal tract (fire and explosion hazard)
- in acute emergency situations
- in the case of severe coagulation disorders
- in the case of severely impaired pulmonary or cardiovascular functions
- Do not use in the area of the heart and central nervous system.

Safe handling and preparation

- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge or experience.
- Read, follow and keep the instructions for use.
- Use the product only in accordance with its intended use, see Intended use.
- Remove the transport packaging and thoroughly clean the new product, either by hand or by a mechanical process, prior to its initial sterilization.
- Store any new or unused products in a dry, clean and safe place.
- Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Replace any damaged components immediately with original spare parts.

Danger of injury from burns and/or explosion from flammable gases!
Application of the HF device according to its intended use can involve spark formation.

- Observe the safety notices in the instructions for use of the HF device.

WARNING

Thermal injuries to patients/users due to insufficient insulation of leads to active accessories

- Adjust the HF device to a setting appropriate to ensure that the maximum peak output voltage does not exceed the accessory voltage rating specified for the product.

WARNING

- Adjust the HF power output to the intended operation. Take into account clinical experience or reference values.
- Select the lowest possible HF power output.
- Keep the product's contact surfaces clean during the operation. Wipe off encrusted tissue residues or bodyfluids, using a moistened swab.

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The rated accessory voltage of the product is 600 V_p.

The accessory voltage rating must exceed or match the maximum peak output voltage at which the product is operated in combination with an appropriate HF device at an appropriate operating mode/setting (see IEC/DINEN60601-2-2).

To avoid HF burns:

- The product's working tip must be in the user's field of vision whenever HF power is activated.
- Prior to activating the HF device, make certain that the working end of the product is not touching any electrically conductive accessories.
- Prior to each use, visually check the product for: damage to, or surface changes on, the insulation.
- Never put down the product on or next to the patient.
- When using accessories for endoscopy or laparoscopy applications, deactivate the automatic switch-on mode of the HF generator.
- Follow the instructions for use of the HF device.

Safe operation

Risk of injury and/or malfunction!

- **Always carry out a function check before using the product.**

WARNING

Validated processing procedure

Note

Adhere to national statutory regulations, international standards and guidelines, and local, clinical hygiene instructions for sterile processing.

Note

Mechanical processing should be preferred over manual cleaning because of the better and more reliable cleaning results of mechanical processing.

Note

Successful processing of this medical product can only be ensured if processing is performed through a validated processing procedure. The user/processor is responsible for the validation.

General notes

To avoid unnecessary, excessive contamination of the complete instrument tray during operations, take care that contaminated instruments are collected separately and not put back into the instrument tray.

Encrusted or fixated residues from surgery can make the cleaning process more difficult or ineffective, and can cause corrosion of stainless steels. Therefore the time interval between application and cleaning should not exceed 2h, and neither potentially fixating pre-cleaning temperatures >45 °C nor any fixating disinfectants (active ingredient:aldehyde, alcohol) be used.

Excessive doses of neutralizers or basic detergents can cause chemical degradation and/or fading and obliteration of inscriptions on stainless steel surfaces, regarding visual reading and machine-readability of the inscriptions. Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfecting and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. To remove such residues, the products must be rinsed sufficiently withfully desalinated water and dried thoroughly.

Only use process chemicals that have been tested and approved (e.g. VAH/DGHM or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations. All process parameters specified by the chemical's manufacturer, such as temperatures, concentrations and exposure times, must be strictly observed.

Failure to do so can result in the following problems:

- Optical deterioration of materials, e.g. fading or discoloration of titanium or aluminum surfaces. For aluminum, pH >8 in the application/process solution can already cause visible surface changes.
- Material damage, e.g. corrosion, cracks, fracturing, premature aging or swelling

- Clean the product immediately after use.

Immersion treatment in a 3 % H₂O₂ solution for approx. 5 minutes is a particularly effective and gentle method to dissolve encrustations from HF instruments. Subsequently, the debris can be removed by hand, with a medium-hard brush and/or in an ultrasonic cleaning bath. After that treatment, proceed through the normal steps of processing. For further detailed advice on hygienically safe and material-/value-preserving reprocessing, see www.a-k-i.org

Preparations at the place of use

- Rinse surfaces inaccessible to visual inspection, e.g. on products with hidden crevices or lumens or products with complex geometries, preferably with distilled water, using e.g. a disposable syringe.
- Remove visible surgical residues as completely as possible, using a lint-free wet wipe.
- Put the dry product into a closed disposal container and have it transferred to cleaning and disinfecting within 2h.

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Cleaning/Disinfecting

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning/disinfecting agents according to manufacturers instructions. The cleaning/disinfecting agents must
 - be approved for e.g. aluminium, plastics and high-grade steel.
 - no attack softeners (e.g. silicone).
- Observe specifications regarding concentration, temperature and exposure time.

CAUTION

- Carry out ultrasound cleaning:
 - as an effective mechanical supplement to manual cleaning/disinfecting.
 - as a pre-cleaning procedure for products with encrusted residues, in preparation for mechanical cleaning/disinfecting.
 - as an integrated mechanical support measure for mechanical cleaning/disinfecting.
 - for additional cleaning of products with residues left after mechanical cleaning/disinfecting.

Manual ultrasonic cleaning and immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Konz. [%]	Water- quality	Chemistry
I	Pre wash	RT (cold)	3	-	T-W	-
II	Disinfecting cleaning	RT (kalt)	15	2	T-W	Borer Chemie; deconex 53 INSTRUMENT, midalkalisch, aldehydfrei, Product data sheet, Borer Chemie: www.borer.ch
III	Intermediate rinsing	RT (cold)	1	-	T-W	-
IV	Final rinsing	RT (cold)	1	-	VE-W	-
V	Drying	RT	-	-	-	-

T-W: Tap water (drinking water)
 VE-W: completely desalinated water (demineralized)
 RT: Room temperatur

Stage I

- Clean the product in an ultrasound cleaning bath (frequency 35 kHz). Make sure that all accessible surfaces are moistened and acoustic shadows are avoided.
- Clean the product with a suitable cleaning brush until all visible residues have been removed from the surface. Brush through all surfaces not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometries, for at least 1 min or until no more residues can be removed, using a suitable cleaning brush.
- Mobilize non-rigid components, e.g. set screws, links, etc. during cleaning.
- After cleaning, thoroughly rinse through these components (at least five times) with the cleaning solution, using a disposable syringe (20ml).
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Stage II

- Rinse the product completely (all accessible surfaces) under running water.
- Allow water to drip off for a sufficient length of time.

Stage III

- Fully immerse the product in the disinfecting solution.
- Rinse lumens at least 5 times at the beginning of the exposure time, using a disposable syringe (20ml) and an appropriate rinsing adapter. Make sure that all accessible surfaces are moistened.

Stage IV

- Rinse the product completely (all accessible surfaces) under running water.
- Rinse lumens at least 5 times, using a disposable syringe (20ml) and an appropriate rinsing adapter.
- Allow water to drip off for a sufficient length of time.

Stage V

- Thoroughly dry the product with lint-free tissue or medical-quality filtered compressed air.

Mechanical cleaning/disinfecting with manual pre-cleaning

Phase	Schritt	T [°C/°F]	t [min]	Wasser- qualität	Chemie
I	Pre- rinse	>25/77	3	T-W	-
II	cleaning	55/131	10	VE-W	Borer Chemie; deconex POWER ZYME Product data sheet, Borer Chemie : www.borer.ch
(III)	(Neutralization, only when cleaning with alkaline cleaning products.	20/68	2	VE-W	- Concentrate, sour: Base: Lemon or phosphoric acid, instructions for use comply with Detergent manufacturer!

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T-W: Tap water (drinking water)
 VE-W: completely desalinated water (demineralized)

Note

The disinfectant must be of tested and approved effectiveness (e.g. DGHM or FDA approval or CE mark according to DIN EN ISO15883).

Note

Use appropriate (neutral, enzymatic and mild-alkaline) detergents for this product with aluminum components.

Note

For thermal disinfection, always use fully desalinated (demineralized) water. Ensure that Ao is >3 000 for the process.

Note

The disinfectant used for processing must be serviced and checked at regular intervals.

Stage I

- Clean the product in an ultrasound cleaning bath (frequency 35kHz). Make sure that all accessible surfaces are moistened and acoustic shadows are avoided.
- Clean the product with a suitable cleaning brush until all visible residues have been removed from the surface.
- Brush through all surfaces not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometries, for at least 1 min or until no more residues can be removed, using a suitable cleaning brush.
- Mobilize non-rigid components, e.g. set screws, links, etc. during cleaning.
- After cleaning, thoroughly rinse through these components (at least five times) with the cleaning solution, using a disposable syringe (20ml).

- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Stage II

- Rinse the product completely (all accessible surfaces) under running water.

Mechanical neutral or mild-alkaline cleaning and thermal disinfection

- Machine type: Single-chamber washer/disinfectant without ultrasound.
- Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).
- Connect components with lumens and channels directly to the rinsing port of the injector carriage.
- Keep working ends open during cleaning.

Inspection, maintenance and checks

Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!

- Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization).

CAUTION

- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the product is: dry, clean, operational and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn or fractured components).
- Dry the product if it is wet or moist.
- Repeat cleaning and disinfecting of products that still show impurities or contamination.
- Check the product for proper functioning.
- Immediately sort out damaged or inoperative products and have them sent to Ulrich Technical Service.
- Check for compatibility with associated products.

Generator connections:

Bipolar forceps from Ulrich Swiss are approved for the use with the following generators:

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COVIDIEN

Type	Item no.	Max. bipolar coagulation power
Force FX™	Force FX™	≤ 70 Watt
ForceTriad™	ForceTriad™	≤ 95 Watt
Force EZ™	Force EZ™	≤ 70 Watt
SurgiStat™	SurgiStat™	≤ 30 Watt

BOWA-electronic GmbH & Co. KG

Type	Item no.	Max. bipolar coagulation power
ARC 400	900-400	≤ 350 Watt
ARC 350	900-351	≤ 350 Watt
ARC 303	900-303	≤ 120 Watt
ARC 250	900-250	≤ 120 Watt
ARC 100	900-100	≤ 100 Watt

Olympus Surgical Technologies Europe

Type	Item no.	Max. bipolar coagulation power
ESG-100, 220-240 V	WB991036	≤ 120 Watt
ESG-100, 100-120 V	WB991046	≤ 120 Watt
ESG-400	WB91051W	≤ 320 Watt

Söring GmbH

Type	Item no.	Max. bipolar coagulation power
MBC 600	MBC 600	≤ 100 Watt
MBC 601	MBC 601	≤ 100 Watt
MBC 601 UAM	MBC 601 UAM	≤ 100 Watt
ARCO 3000	ARCO 3000	≤ 100 Watt
MBC 200	MBC 200	≤ 140 Watt
BCC 140	BCC 140	≤ 140 Watt

LAMIDEY NOURY MEDICAL

Type	Item no.	Max. bipolar coagulation power
Optima 4	Optima 4	≤ 370 Watt
Optima 3	Optima 3	≤ 370 Watt
Optima 2	Optima 2	≤ 200 Watt
MC 2	MC 2	≤ 70 Watt
MC 3	MC 3	≤ 100 Watt
MC 4	MC 4	≤ 100 Watt

Integra LifeSciences Corporation

Type	Item no.	Max. bipolar coagulation power
Elektrotom® 621	Elektrotom® 621	≤ 95 Watt
Elektrotom® 630	Elektrotom® 630	≤ 95 Watt

Gebrüder Martin GmbH & Co. KG

Type	Item no.	Max. bipolar coagulation power
Minicutter	80-008-03-04	≤ 70 Watt
ME 102	80-010-02-04	≤ 80 Watt
maxium®, m-Version	80-042-00-04	≤ 320 Watt
maxium®, i-Version	80-042-02-04	≤ 320 Watt
maxium®, e-Version	80-042-04-04	≤ 320 Watt
ME MB 3, m-Version, 220-240 V	80-040-11-04	≤ 100 Watt
ME MB 3, m-Version, 100-127 V	80-040-11-10	≤ 100 Watt
ME MB 3, i-Version, 100-127 V	80-040-12-04	≤ 100 Watt
ME MB 3, i-Version, 100-127 V	80-040-12-10	≤ 100 Watt
ME 411	80-041-01	≤ 80 Watt

Aesculap AG / B. Braun Melsungen AG

Type	Item no.	Max. bipolar coagulation power
GN300	GN300	≤ 80 Watt
GN640	GN640	≤ 100 Watt
GN060	GN060	≤ 50 Watt

Suitable connection cables:

Ulrich Swiss, bipolar forceps have been approved for use with the following bipolar connection cables for the above generators:

Cable connection, generator
ERBE Elektromedizin GmbH, KARL STORZ GmbH & Co. KG
ERBE Elektromedizin GmbH – ICC-Serie, International
Aesculap AG / B. Braun Melsungen AG, Gebrüder Martin GmbH & Co. KG, Sutter Medizintechnik GmbH
COVIDIEN, LAMIDEY NOURY MEDICAL, BOWA-electronic GmbH & Co. KG, ERBE Elektromedizin GmbH - serie VIO, Olym-pus Surgical Technologies Europe, Söring GmbH
Cable connection, forceps
European flat plug
2-pin plug (US male connector)

These cables may be connected only to the bipolar output of electrosurgical devices. Please follow the generator instructions on connecting these cables before using the bipolar forceps. HF generators pose a number of risks, such as incorrect operation, unintended high-frequency burns, ignition of flammable liquids and gases (risk of explosion).

Packaging

Appropriately protect products with fine working tips.

- Sort the product into its appropriate storage device or put it on a suitable tray. Make certain that all cutting edges are protected.
- Package trays appropriately for the sterilization process.
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage (DIN EN ISO 11607).

Sterilization method and parameters

- Make certain that all external and internal surfaces will be exposed to the sterilizing agent (e.g. by opening all valves and faucets).
- Validated sterilization process
 - Steam sterilization through fractionated vacuum process
 - Steam sterilizer according to DIN EN 285, validated according to DIN EN ISO17665
 - Sterilization through fractionated vacuum process at 134 °C/holding time 5 min
- When sterilizing several products at the same time in one steam sterilizer: Make certain that the maximum load capacity of the steam sterilizer, as specified by the manufacturer, is not exceeded.

Storage

- Store sterile products in germ-proof packaging under dust protection in a dry, dark and temperature-controlled room.

Technical Service

Risk of injury and / or malfunction!

- **Do not modify the product**

Warning

For service and repair, contact Ulrich AG, St. Gallen. Modifications to medical equipment may result in loss of Guarantee / warranty claims as well as loss of medical approvals.

Service Address

Ulrich AG
Mövenstrasse 12
CH – 9015 St. Gallen

Note: This instruction for use was written in German and translated into English, French and Italian.

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